Helicobacter pylori AMA Rapid Urease Test Expert

AMA RUT Expert with AMA RUT Reader

INSTRUCTIONS FOR USE

English.....	4
Español.....	25
Français.....	51
Deutsch.....	75
Suomalainen.....	100
Svenska.....	121
Română.....	145
Ελληνικά.....	169
العربية.....	223

4

Thank you for choosing AMA RUT Expert, the rapid urease test for semi-quantitative detection of *H.pylori*. Please follow the instructions carefully, and the test-system will help to diagnose your patients better and save your time.

INTENDED USE AND PRINCIPLE OF OPERATION

The test-system consists of a single use test AMA RUT Expert (hereinafter referred to as slide) and the AMA RUT Reader. The intended use is specific rapid *H.pylori* detection by establishing the presence of urease activity in a biopsy specimen taken from either adult or child patients during the gastroscopy procedure.

Biomaterial test samples could be:


- a biopsy specimen taken from any part of the stomach;
- a biopsy specimen taken from the duodenal cap.

One to three specimens can be tested at the same time. The size of a biopsy specimen should be no less than 2 mm at any dimension.


The principle of operation is based on the colour change of the indicator after the biopsy specimen(s) placed on the reactive element of the slide. In the event of urease activity in the biopsy specimen(s), a red or magenta spot appears on the slide. The Reader takes multiple readings to check if the spot appeared. The higher the urease activity is, the faster the detection time. The Reader stops as soon as it detects a spot or proceeds up to the maximum exposure time (14 minutes) to confirm the absence of the spot and thus the negative result of the test. The result is indicated on the display and saved.

DESIGN OF THE TEST-SYSTEM

AMA RUT Expert is a rectangular-shaped polymer slide with the reactive element hermetically sealed by a protective cover.


The Reader is a cube-shaped electronic device with the LCD display and the operating button on the top side, the battery compartment and the USB cable socket are on the lateral side. Slide holder with the slot situated on the bottom side.

6

The Reader can be powered by batteries (3 pcs., 3V CR2032) or from the mains via an **AC/DC USB adapter*** (optional, output 5 V, 500 mA) and a specific **USB cable*** (optional).

PREPARATION BEFORE THE TEST

- Put on the gloves.
- Open the protective cover gently. Do not tear it off.
- Put the slide on a flat surface.

THE PROCEDURE

1. Use dry clean forceps or the biopsy pick. Place the biopsy specimen(s) on the reactive element of the slide. There is a place for up to three specimens on the reactive element.

**The biopsy specimen(s)
should be placed directly on
the reactive element of the
slide and should not extend
beyond its borders.**

CD with the software and driver*
(optional)


Control sample* (optional)

MATERIALS REQUIRED, BUT NOT PROVIDED

- Forceps
- Powder-free gloves

2. When the last biopsy is placed, reseal the cover. Squeezing the seal not recommended, but accurate resealing is important to prevent the biopsy specimen(s) from drying up.


3. Turn the slide upside down so that the yellow reactive element positions above.

4. Insert the slide into the slot on the bottom side of the Reader as far as it can go.

The slide should be inserted with the yellow reactive element and the arrows upside. Make sure the period between the biopsy placement onto the slide and the first reading is as short as possible. Should the period exceed 1 minute, the test result is interpreted as qualitative only, please see *Limitations* on page 20.

5. Turn the Reader on by pressing the button. When the display check is completed, the display shows "WAIT" following; a "QC99" message appears. Alternatively an audible alarm may beep three times and the display shows some default date and time values (for example, "2016.03.21 11:17"), it means the correct date and time settings were lost. However, the test can be continued, but the result will be saved under the default values of time and date. Before the next test, please follow the instructions for date and time setting on page 13.

QC99

"QC" stands for "quality check", the procedure that must be fulfilled once in every 1000 tests to check the correct reading. When the number of the tests comes close to the check, the message turns

“QC98”, “QC97” etc. until it reaches “QC00”. At this point, no more readings can be taken until the quality check procedure is passed successfully.

Then, if previous readings have been taken, an audible alarm beeps and the display shows last saved result as a scrolling message. Press the button to continue, the audible alarm beeps and the display shows “ON”.

6. To start the test, press the button again. The display shows “TEST”.

7. Confirm the start of the test by pressing the button. The Reader first checks if the slide is properly inserted into the slot and if yes, the first reading is taken and the display shows “RUN”.

If instead of “RUN” the display shows “QC-TEST-NEEDED”, please refer to *Performing the QC procedure on page 16*.

If the display shows “<---”, it means the Reader cannot find the slide in the slot. Make sure the slide is inserted as far as it can go and the result evaluation side up. Press the button to restart the procedure from step 6.

After the first reading, the Reader either stops the procedure with the positive result, or continues with a countdown showing time until the next reading. After the counter reaches zero, the next reading is taken, and the display shows “RUN” again. The countdowns and the readings repeat automatically until the slide is interpreted as positive or until total time reaches 14 minutes. The display then shows the result as a scrolling message.

ON

TEST

RUN

<---

00:25

Period of exposure since the biopsy placement	Result on the screen	Urease activity appropriate for the number of bacterial <i>H.pylori</i> cells*
0 and 1 st minute	HP +++	High, as in ~10 ⁷ CFU and more
3 rd to 4 th minute	HP ++	Moderate, as in ~10 ⁶ CFU
9 th to 14 th minute	HP +	Low, as in ~10 ⁵ CFU

*The approximate values are for indicative purposes only, these are based on the microbiological research, see *References* on page 23. AMA RUT Expert does not in any way measure the number of the *H.pylori* cells.

Inserting the slide into the Reader with delay from the moment of biopsy placement increases the exposure time, thus the low or moderate urease activity cases may indicate as the “HP+++” results.

The “1”, “2”, “3” in the scrolling message of the positive indicated result specify the biopsy placement area(s) where the spot has appeared the fastest, for example “HP++ 1,2”.

The absence of the spot on the 14th minute means the absence of the urease activity. The “HP-” on the screen indicates the negative test result.

The “Temp<17” in the scrolling message means the room temperature was below 17 °C. In such cases, the speed of the enzymatic reaction is slowed down and the low urease activity might not be detected. It is recommended to repeat the reading of the same slide in normal room temperature, starting from step 3 on page 7, and interpret the result as qualitative only.

After completing the test, the biopsy specimen is applicable for further examinations, such as PCR, histology or culture.

There is no button to switch the device off. The Reader continues beeping for 3 minutes, and then shuts down automatically if no actions are made. The results are saved and shown as a scrolling message at next switch-on.

TRANSFERRING THE TEST RESULTS

The Reader memory can store 86 test results. When the memory is full, old data is replaced with the new.

To avoid data loss, transferring the test results from the Reader to the PC is recommended weekly.

To transfer the test results, the Reader has to be successfully connected to the software, please see *Installing the software* on page 14.

1. Click the *Get data from Reader* button. The software will receive the entries from the memory of the Reader and display them in the main window. The window is divided into two parts. The first part, on the left, lists all the test results. The list contains the *Test ID* (a unique number added by the Reader to each test result), *Test name*, *Lot number*, *Manufacturer*, *Date/Time* and the *User comments* (can be added manually for every entry listed). It is possible to sort the entries in alphabetical or value dependant matter. This can be done with every column. By clicking on one of the entries, the actual result values for it will be presented in the list on the right side of the window. This second list will

show the actual *Result* and *Class* as well as some additional technical data, please see an example on the right.

The first row of the *Result* indicates the period of exposure since the biopsy placement, in this example it is "2,3" minutes (135 seconds). The fifth row of the *Class* indicates the actual test result, in this example "HP++". The sixth row of the *Class* indicates the area, where the spot appeared, in this case the one in the middle, "2". In case the spots are detected in all the three areas, the "1,2,3" value is indicated in the seventh row of the *Class*.

The next example shows a negative test result:

The period of exposure is "14,0" minutes (the first row of the *Result*), the test result is "HP-" (the fifth row of the *Class*) and the room temperature is "24,5" °C (the second row of the *Result*).

Result	Result unit	Class
2,3		
500,0		
-18,0		
4,0		
		HP++
		2
		1,5

Result	Result unit	Class
14,0		
24,5		
-68,0		
14,0		
		HP-
		100,0

 **The other values presented in the list on the right side of the window are technical and not intended for customer use. Please contact the manufacturer for more details, if additional explanations are needed.**

2. To save the transferred test results click on *Menu* and *Save data to file* in the menu bar. The saved file with ".tsv" extension (tab separated values) can be opened by AMA RUT Reader software and by other applications, for example Microsoft Excel.

When more test results have been completed and saved to the same file, a pop-up window appears, saying "Do you want to add the data to the selected file? If not, the file will be overwritten with the

new data." Press "yes" to update the existing file with the new data. To have the database indicated at next launch of *AMA RUT Expert* software, click on *Settings* and *Load last used file on start*.

3. After transferring the test results is finished, press the *Disconnect* button and then unplug the Reader from the PC.


CHANGING THE BATTERIES

The status of the batteries is automatically monitored. Fresh batteries of well-known brands allow about 70 tests to be performed. Replace the batteries when the battery symbol starts to blink after a test. Three fresh batteries CR2032 are required. Remove from the packaging and insert into the battery compartment of the device.

i We recommend changing the batteries when the Reader is connected to the mains via the USB cable. In this case, the date and time settings are not lost while the batteries are out. Without the cable connection, the time and date settings are lost if the batteries don't power the Reader for 1 minute or longer.

Turn the battery cover with a smooth-edged coin counter-clockwise until it stops. Turn the Reader slightly down, removing the battery cover. Discard the used batteries, they cannot be recharged and have to be properly disposed. Place the three fresh batteries with correct polarity ("+" side up) into the tray one by one. Thereafter, the battery cover has to be slightly turned with a coin and rotated clockwise until it stops.

While inserting the batteries body grease attached to fingers should not come in contact with the batteries. Contamination can lead to a more rapid discharge of the batteries. Therefore gloves or

plastic tweezers are recommended. In case the Reader does not start with the fresh batteries please check the polarity and clean the batteries by using a dry cloth.

SETTING DATE AND TIME

The date and time settings of the Reader set at manufacture and are stored as long as the batteries power the Reader. It is important to keep the settings correct because the date and time help to distinguish between the results stored in the internal memory. The settings are lost if the Reader is not powered for 1 minute or longer. In this case, it is necessary to set the date and time.

The easiest way to set the date and time is to connect the Reader to a PC with the USB cable when the AMA RUT Expert software is running (please follow the instructions on the software installation on page 14). The correct settings will be applied automatically as soon as the connection is established.

If the software is not available, we recommend setting the date and time when the Reader is powered from the mains via the USB cable and the adapter. When powered from the mains, the Reader does not shut down automatically and more time for the setting is available.

If neither the software nor the cable is available, please do the following within 2 minutes:

1. Bring the device to position “ON” by pressing the button. Press the button longer than 1 second, but no longer than 5 seconds, and then immediately one more time. Year, date and time will appear on the display.

Make sure the display doesn't show “TEST”. In this case, press the button few times until the display shows “ON” and try again.

14

2. Press the button longer than 5 seconds, a flashing display will show the first time specification, year.
3. By repeated short (<1 second) pressing the button, change the displayed value. When the wanted value is reached (e.g. year) press the button longer (>5 seconds), the appointed value will be stored and the next value will be presented.
4. Repeat these steps to set to year, month, day, hour and minute. After setting the date and time accordingly the display will show "OK".
5. Press the button one more time, the display will show "ON" and the Reader is now ready for use.

If the Reader remains without power for 1 minute or longer, the display will show year, date and time with default values at first switch-on. In this case, start the setting from the step 2.


INSTALLING THE SOFTWARE

To install the software, launch the setup executable from the provided CD. Please follow the instructions of the installer. After the installation is finished, it is recommended to keep the checkmark for FTDI driver installation and to proceed with it. This is necessary to have the connection between the Reader and the PC established.

The system requirements to install the software are:


- Operating system: Windows XP, Vista, 7 (SP1), 8, 10
- Processor: Pentium III or higher, AMD Athlon XP or higher
- Main memory: min. 1 GB
- Hard disk: min. 1 GB free
- Screen resolution: min. 1336 x 768

After the installation has been successful, press the button of the Reader until its display shows "ON", connect the Reader to the PC via the provided USB cable, then launch the *AMA RUT Expert* software.

When the Reader is connected to the PC via the USB cable, it is powered by the USB port and doesn't discharge the batteries. There is no auto-shutdown until the Reader is unplugged.

As soon as the Reader is ready for operation the software will automatically try to connect to the Reader. This automated process operates preconfigured after installing the software and can be switched off later. In such case, it is possible to connect the Reader manually by pressing the *Connect* button. If the Reader is connected successfully, the button turns to *Disconnect*.

The software activates automatically. To check that, click on *Info* and *About AMA RUT Expert*. The window will say that the software license is valid.

To set the interface language that suits you best use the *Language* in the menu bar on top of the software window. The Reader is now connected and ready to work with the software.

PERFORMING THE QC PROCEDURE


For quality control purposes, the correct reading of the device must be checked through the quality check procedure (QC). The Control samples are being read and the results are being compared to the nominal values that are stored in the special QC configuration file from the provided CD.

The procedure can run at any time as a facultative quality control measure, but when the number of the tests performed comes to 1000, the quality check becomes mandatory. The Reader is locked, and as soon as the procedure has been successfully passed, the counter clears to zero and another 1000 tests can be performed.

Before the first quality check it is necessary to specify the location of the QC configuration file. Copy the *QC configuration.bin* file from the provided CD to the hardware of your PC. In the menu bar on top of the software window click on *Menu* and *Select QC file*, a window will pop up where you have to state the location of the *QC configuration.bin* file on the hardware.

To perform the quality check, make sure the Reader is connected and the software is running as described in *Installing the software* on page 14 and do the following:

1. Press the *Quality Check* button. A pop-up window will appear, saying "Insert the Control sample into the AMA RUT Reader".
2. Insert the Control sample into the slot on the bottom side of the Reader as far as it can go and press the *OK* button in the pop-up window.

After the reading of the Control sample is done, the QC result ("Pass" or "Fail") is shown both on the display of the Reader and in the "Class" column on the right side of the window.

Result name	Result	Result unit	Class	Result name	Result	Result unit	Class
QC			Pass	QC			Fail

If the QC result is "Pass", the Reader can be used further. If the QC result is "Fail", please try to do the following:

- 1) If the "Result" column is shown as "<----", it means the Reader cannot find the Control sample in the slot. Make sure the Control sample inserts as far as it can go. Then press the *Quality Check* button to restart the procedure.
- 2) Make sure the Control sample provided by the manufacturer is used for the QC procedure. Please avoid the shorter slide for diagnostics (AMA RUT Expert) to be used for the QC procedure.
- 3) Detach the slide holder and clean the glass on the bottom of the Reader (please see *Maintenance and cleaning of the Reader* on page 21).
- 4) Make sure the Control sample is not dirty.
- 5) If none of these helps, contact the Manufacturer and be ready to provide with the saved QC result(s).

Validity	Result name	Result
Valid	QC	
Invalid		<----

3. To save the QC test result, click on *Menu* and *Save data to file* in the menu bar. The saved file with ".tsv" extension (tab separated values) can be opened by AMA RUT Reader software and also by other applications, for example Microsoft Excel.

We strongly recommend saving the results of every QC procedure. If the QC result is not saved before closing the software, it will be lost.

OTHER OPTIONS OF THE SOFTWARE

To have all options selectable, The Reader needs to be turned on, connected to the software, and the test results from the Reader memory should be transferred to the software.

When clicking on *Menu* the following options are available:

- *Delete results list*: By choosing this option, all entries in the current list delete.
- *Save data to file*: By choosing this option, a window will open which allows picking a place to save the reading results as a *.tsv file (tabular separated values). Such a file can be directly opened with Excel, all information will be separated into different columns.
- *Load data from file*: By choosing this option, a previously saved results list can be reloaded. The list has to be a *.tsv file.
- *Format AMA RUT Reader database*: By choosing this option, the Reader's database/memory will be formatted. All data will be lost.

When clicking on *Settings* the following options are available:

- *Set image position*: When doing a quality check, this option determines where to place the picture of the 2D volume diagram in conjunction with the results list. Presented will be the volume diagram from the actual reading.

Possible image positions are the following: on the left side of the window, on the right side of the window, full width of the window, do not display the image.

- *Data management*: This option allows choosing how the software will proceed with new data either just read or downloaded, if there is already data in the list.

The following options can be selected:


- Overwrite existing data: by choosing this option, all entries in the list will be deleted when adding new data either by doing a quality check or by adding data from another reader.
- Add new data to existing: by choosing this option, all entries in the list will remain when adding new data either by doing a quality check or by adding data from another reader.
- Always ask: the user will always be asked how to manage new data. The user can choose between adding the data and overwriting the existing data.
- *Auto connect*: by choosing this option, the user can decide whether the software should automatically connect to the Reader or if the connection should be done manually.
- *Shutdown AMA RUT Reader on disconnect*: by choosing this option, the Reader switches off when being disconnected from the software.
- *Load last used file on start*: by choosing this option, the recently saved results list restores when restarting the software.

WARNINGS AND PRECAUTIONS

For In Vitro Diagnostic Use. For use by trained personnel only.

Handle biopsy specimens as potentially biohazardous material. All biopsy specimens should be regarded as potentially contaminated and treated as if they were infectious. Please refer to the local or national regulations – see References on page 23. Always use protective gloves when handling patient samples. Read all instructions prior to performing the test. Do not use the slide beyond the expiry date. Discard the used slides to biohazardous waste according to the local and national regulations.

All Readers produced, and checked before shipping under strict quality control measures in order to guarantee a high degree of quality. Do not attempt to re-configure or open the Reader (except detaching the slide holder). Always ensure that the slide holder is attached to the Reader correctly. A wrong or improper attachment can lead to incorrect results. Do not use the Reader in direct sunlight or exposed to bright light while the readings are being taken. The Reader is designed for use on a straight and horizontal surface area. The Reader can be operated at temperature range from 17 to 35 °C and humidity range between 20% and 85%. Ensure that the Reader reaches operating temperature before use. Protect the Reader from any liquids. Any liquid entering the Reader enclosure can damage the Reader permanently.

LIMITATIONS

False negative results may occur if:

- *H.pylori* inhibiting antibiotics have been taken 4-6 weeks prior to the examination;
- Acid inhibiting drugs (PPI or H₂ blockers), antisecretory agents, anti-inflammation agents, analgetics and bismuth drugs have been taken 14 days prior to the examination;
- Low ambient temperature; see *The procedure* on page 6.

Semi-quantitative results of AMA RUT Expert are based on kinetic evaluation of the urease activity. For this, it is important to fix the period from the moment the biopsy specimen(s) carefully placed on the reactive element until the possible spot appearance. The time fixation is carried out by the Reader automatically, but it is important to start the reading no longer than 1 minute after the biopsy specimen(s) have been placed on the reactive element of the slide. The semi-quantitative result is valid only in this case. Otherwise, if the delay is longer than 1 minute, the test results should be

interpreted as qualitative only. As with any diagnostic procedure, the test results indicated on the screen are interpreted in the light of the patient's clinical presentation and any other information available to the physician.

MAINTENANCE AND CLEANING OF THE READER

- The Reader is maintenance free.
- Before each test, the glass on the bottom (covered by the slide holder) should be checked for impurities. For cleaning a commercial cloth together with a commercial cleaning fluid, e.g. for glasses products, is recommended.
- Due to the possibility of contact with potentially infectious material while using the Reader, disinfection is required to prevent infections or other contaminations. For disinfection every surface of the Reader as well as the slide holder needs cleaning with a suitable disinfection solution. That solution approved for medical devices, should not affect the device housing material. The disinfection spray Mikrozid® AF Liquid for instance is a suitable solution. Comparable products can be used as well.
- For cleaning the display use dry cloth. Never apply fluid onto the display!

STORAGE, STABILITY AND TRANSPORTATION TERMS

Store the slide and the Control sample

- in the manufacturer's packaging;
- in a dark, dry place with the temperature from +15 °C to +50 °C;
- in a place protected from mechanical actions (friction, pressure, strokes);
- keep away from the ammonia vapor, moisture and direct sunlight.

When stored at this temperature the slide and the Control sample are stable for 24 months. Expiration date of the slide indicated on the package label. Expiration date of the Control sample is indicated right on it. Transport by any kind of transportation with the temperature from -50 °C to +60 °C, sealed. The transport period must not exceed 1 month.

Store the Reader

- in the manufacturer's packaging;
- in a dark, dry place with the temperature from -30 °C to +80 °C and humidity from 20% to 85%.

When stored at this conditions the Reader's lifetime is 5 years. Date of manufacture indicated on the Reader's label.

WARRANTY

The manufacturer shall remedy all defects discovered in any Product (the "Defective Product") that result from unsuitable materials or negligent workmanship and which prevent the mechanical functioning or intended use of the Products including, but not limited to, the functions specified in the manufacturer's specifications for the Products.

Any warranty will, however, be deemed as void if fault is found to have been caused by maltreatment, misuse, accidental damage, incorrect storage or use of the product for operations outside their specified limitation or outside their specifications, contrary to the instructions given in the instruction manual.

The period of the AMA RUT Expert warranty is 24 months from the date of manufacture. The period of the Reader's warranty is 24 months from the date the relevant Product was shipped by the manufacturer.

ORDERING INFORMATION

- AMA-Med Oy, Sammonkatu 12, 50130, Mikkeli, Finland, Tel: +358 45 164 44 04, E-mail: expert@amarut.com
- Manufacturer: Association of Medicine and Analytics Company Limited, 17 line of Vasilievsky Island, 4-6, 199034, St-Petersburg, Russia, Tel: (007) 812 321-7501, E-mail: ama@sp.ru

REFERENCES

1. St. Petersburg Pasteur Institute of Epidemiology and Microbiology, St. Petersburg, Russia. Protocol of laboratory evaluation of AMA RUT Expert with the AMA RUT Reader for the detection of Helicobacter pylori urease with respect to the speed of reaction and detection limits, 2016.
2. St. Petersburg State Budgetary Healthcare Institution "City Mariinsky hospital", St. Petersburg, Russia. Protocol of clinical evaluation of AMA RUT Expert with the AMA RUT Reader for the detection of Helicobacter pylori infection in biopsy specimens, 2016.
3. U.S. department of Health and Human SERVICES (Bethesda, MD., USA) publication Biosafety in Microbiological and Biochemical Laboratories, 1999, 4th ed. (CDC/NIH) and No. (CDC) 88-8395 on reports of laboratory safety procedures on different diseases.

EXPLANATION OF THE SYMBOLS, USED IN LABELS

Manufacturer

Use by

Date of manufacture

Do not reuse

The Reader is to be disposed in accordance with the applicable country-specific regulations

In vitro diagnostic medical device

European Conformity

Authorised Representative in the European Community

Batch code


Contents sufficient for N tests

Consult Instructions for use

Temperature limitation

Serial number

AMA is the internationally registered trademark of AMA Co Ltd. Windows, Microsoft, Excel, Athlon, FTDI, Pentium, AMD, Vista, Mikrozid and etc. are registered trademarks that belong to their respective owners.

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LÄHTEET

3) St. Petersburg Pasteur Institute of Epidemiology and Microbiology, St. Petersburg, Russia. Protocol of laboratory evaluation of AMA RUT Expert with the AMA RUT Reader for the detection of *Helicobacter pylori* urease with respect to the speed of reaction and detection limits, 2016.

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4) St. Petersburg State Budgetary Healthcare Institution "City Mariinsky hospital", St. Petersburg, Russia. Protocol of clinical evaluation of AMA RUT Expert with the AMA RUT Reader for the detection of *Helicobacter pylori* infection in biopsy specimens, 2016.

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
REFERINȚE

1. St. Petersburg Pasteur Institute of Epidemiology and Microbiology, St. Petersburg, Russia. Protocol of laboratory evaluation of AMA RUT Expert with the AMA RUT Reader for the detection of Helicobacter pylori urease with respect to the speed of reaction and detection limits, 2016.
2. St. Petersburg State Budgetary Healthcare Institution "City Mariinsky hospital", St. Petersburg, Russia. Protocol of clinical evaluation of AMA RUT Expert with the AMA RUT Reader for the detection of Helicobacter pylori infection in biopsy specimens, 2016.

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