Helicobacter pylori

PRUEBA RÁPIDA DE UREASE Expert AMA (AMA RAPID UREASE TEST Expert) (AMA RUT Expert)

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14. NORMATIVAS

1. Mobley H.L.T., Mendz G.L., Hazell S.L., editors. *Helicobacter pylori: Physiology and Genetics*. Washington (DC). ASM Press. 2001.
2. Ji-Hyun Seo, Ji-Sook Park, Kwang-Ho Rhee, and Hee-Shang Youn. Diagnosis of *Helicobacter pylori* Infection in Children and Adolescents in Korea. *Pediatr Gastroenterol Hepatol Nutr*. 2018 Oct; 21(4): 219–233.
3. Lash J.G., Genta R.M. Adherence to the Sydney System guidelines increases the detection of *Helicobacter gastritis* and intestinal metaplasia in 400 738 sets of gastric biopsies. *Alimentary Pharmacology & Therapeutics*. 2013. Aug; 38(4): 424-431.
4. Mollenkopf C, Steininger H, Weineck G, Meyer M. Gastritis: immunohistochemical detection of specific and nonspecific immune response to *Helicobacter pylori*. *Z Gastroenterol*. 1990 Jul; 28(7): 327-34.
5. Federal State Budgetary Institution "Federal Scientific and Clinical Center for Specialized Types of Medical Aid and Medical Technologies of the Federal Medical and Biological Agency", Clinical Protocol Evaluation Protocol, Russia, 2019.
6. Saint-Petersburg Pasteur Institute of Epidemiology and Microbiology, St. Petersburg, Russia. Study Protocol No. 6/2, 2017.
7. Lash J.G., Genta R.M. Adherence to the Sydney System guidelines increases the detection of *Helicobacter gastritis* and intestinal metaplasia in 400 738 sets of gastric biopsies. *Alimentary Pharmacology & Therapeutics*: 2013 Aug; 38(4): 424-431.
8. Bode G., Mauch F., Malfertheiner P. The coccoid forms of *Helicobacter pylori* – criteria for their viability. *Epidemiol. Infect.* 1993; 111: 483-490.
9. Cellini L., Allocati N., Angellutti D., et al. Coccoid *Helicobacter pylori* not culturable in vitro revers in mice. *Microbiol. Immunol.* 1994; 38: 834-850.
10. Takahiro Uotani, David Y. Graham. Diagnosis of *Helicobacter pylori* using the rapid urease test. *Ann Transl Med*. 2015. Jan; 3(1): 9.

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EN

Helicobacter pylori AMA RAPID UREASE TEST Expert (AMA RUT Expert)

Express-test for the Detection of *Helicobacter pylori* Urease Activity in Biopsy Specimen. For in vitro diagnostic use. Store at 4-42 °C Upon Receipt.


INSTRUCTION FOR USE

Authorized Representative in Europe "Emergo Europe", Prinsessegracht 20, 2514 AP The Hague, The Netherlands	Manufacturer "Association of Medicine and Analytics Company Limited", 17 line of Vasilevsky Island, 4-6, 199034, Saint-Petersburg, Russia	

EDITION 05, December, 2020

1. INTENDED USE AND PRINCIPLE OF OPERATION

The intended use is specific rapid *Helicobacter pylori* detection by establishing the presence of urease activity in a biopsy specimen taken from either adult or child patients during the gastroscopy procedure. *Helicobacter pylori* produces the enzyme urease, which breaks down urea to ammonia. The detection of the enzymatic activity is based on the following biochemical reaction [1]:

The principle of operation of AMA RUT Expert is based on the color change of the indicator after the biopsy specimen(s) placed on the reactive element of the test. In the event of urease activity in the biopsy specimen(s), a red or magenta spot appears on the back side of the test. The sensitivity threshold is 2.6×10^4 CFU *Helicobacter pylori*. Biomaterial test samples could be:

- A biopsy specimen taken from any part of the stomach;
- A biopsy specimen taken from the duodenal cap [2, 3, 4].

One, two or three specimens of the one patient can be tested at the same time. The size of a biopsy specimen should be no less than 2 mm at any dimension.

2. DIAGNOSTIC CHARACTERISTICS

For evaluation of clinical performance, the comparison of AMA RUT Expert results was done with the results of histological and/or PCR analysis of biopsies from 200 patients (2 biopsies for each patient). In all the cases one biopsy specimen was placed on the AMA RUT Expert first. Indicative effect was observed after 14 minutes, and then biopsies specimens were sent to histological and/or PCR analysis. Sensitivity was defined as a ratio of the amount of positive results of AMA RUT Expert towards the total amount of HP-positive patients diagnosed with the reference methods. Similarly, specificity was defined as a ratio of negative results according to AMA RUT Expert towards the total amount of HP-negative patients. In that regard, clinical performance of the AMA RUT Expert is the following: sensitivity is 98 %, specificity is 100 %, accuracy is 99 %, negative predictive value is 98 %, positive predictive value is 100 % (Table 1) [5].

Table 1. Diagnostic characteristics

Diagnostic characteristic		True positive		Total
		+	-	
AMA RUT	+	103	0	103
	-	2	95	97
Total		105	95	200
Se, %		98 %		
Sp, %		100 %		
Accuracy, %		99 %		
NPV, %		98 %		
PPV, %		100 %		


Additionally, the cross-sensitivity towards non-target urease-positive bacteria on test results was investigated. The study was conducted using cultures of *Staphylococcus aureus*, *Proteus mirabilis*, *Klebsiella pneumonia* (concentration 10^9 CFU/ml) and the control culture *Helicobacter pylori* (10^4 CFU/ml). The result of the study showed the negative reaction on the non-target bacteria and weak positive on the control culture after 14 minutes [6].

Table 2. Influence of urease-positive bacteria on test results

Control	Approximate number of microorganisms	Result
<i>Helicobacter pylori</i>	10^4	Positive
<i>Proteus mirabilis</i>	10^9	Negative
<i>Klebsiella pneumonia</i>	10^9	Negative
<i>Staphylococcus aureus</i>	10^9	Negative

3. DESIGN OF THE TEST

The test is a rectangular-shaped polymer basement with the reactive element hermetically sealed by a protective cover.

4. DELIVERY SET

Delivery set should satisfy the requirements in the Table 3.

Table 3.

Product name	Quantity, pcs.
<i>Helicobacter pylori</i> AMA RAPID UREASE TEST Expert (AMA RUT Expert)	50 / 48 in blocks / block of 16 / block of 8
Biopsy pick (optional)	50 / 48 / 16 / 8
Instruction for use	1

The block consists of 8 or 16 pieces of the tests, each test can be separated and functions independently. The tests in the block are separated by perforation, allowing to tear a test off and work with it.

5. WARNINGS AND PRECAUTIONS

For in vitro diagnostic use. For use by trained personnel only. After the test is done, the biopsy specimen is applicable for further histology, PCR or culture detections.

Caution: Handle biopsy specimens as potentially biohazardous material.

All biopsy specimens should be regarded as potentially contaminated and treated as if they were infectious. Please refer to the local or national regulations. Always use protective gloves when handling patient samples. Read all instructions prior to performing the test. Do not use the test beyond the expiry date. Discard the used devices to biohazardous waste according to the local and national regulations.

6. MATERIALS REQUIRED, BUT NOT PROVIDED

- Forceps
- Timer
- Powder-free gloves

7. LIMITATIONS

Not recommended:

- testing at an ambient temperature below 17 °C and above 35 °C;
- eating legumes a day before testing;
- drinking carbonated beverages 2-3 hours before testing.

All these factors can affect the reliability of test results, in addition, false negative results may occur if:

- the concentration of in the sample taken during the biopsy is below the sensitivity threshold of the test;
- the sample does not contain *Helicobacter pylori* due to the uneven distribution of bacteria in the gastric mucosa (in 1-5 % of patients, the bacteria is present in the body/corner of the stomach, and not in the antral part, and conversely) [7];
- if the sample contains bacillary forms of *Helicobacter pylori* cells below 10 % of the total number of cells [8, 9];
- the biopsy stick or forceps contaminated before testing;

- taking antibiotics that inhibit *Helicobacter pylori* 4-6 weeks before the test;
- PPIs or H2 receptor antagonists, analgetics, antiseptic agents, anti-inflammation agents, bismuth compounds have an anti-*Helicobacter pylori* activity and decrease the load of *Helicobacter pylori* leading to false negative results, an interval of 2 weeks is required [10].

The test results must be interpreted by the physician based on the patient's clinical presentation. If there is a discrepancy with other diagnostic parameters it is recommended to make additional tests using other methods.

8. PREPARATIONS BEFORE THE TEST

- Put on the gloves
- Open the package of the test


9. TEST PROCEDURE

If you use AMA RUT Expert in blocks, separate one test along the perforation lines and put it on a flat surface.

Open the protective cover, but do not pull hard not to tear it off.

Place one, two or three biopsy specimen(s) on the white reactive element, but not beyond.

After the last biopsy is placed, reseal the cover. Gently smooth the cover with your finger around the perimeter of the test.

If you use AMA RUT Expert for a visual evaluation, check the results from the back of the test after 14 minutes.

10. EVALUATION OF TEST RESULTS

The presence of the red or magenta spot on the reactive element of the test indicates urease activity in the biopsy specimen(s). The greater the urease activity is the lower the detection time will be. If the reactive element remains yellow, then the urease activity in the biopsy specimen(s) is absent. Examples of evaluating test results are shown in Table 4.

Table 4. Matching the observed colors of the reactive element to the field numbers

Colors of the reactive element					Results
1	2	3	4	5	
Yellow	Yellow	Yellow	Yellow	Yellow	Absence of urease activity HP-
Yellow	Yellow	Yellow	Yellow	Red	Presence of urease activity HP+
6	7	8	9	10	

11. STORAGE, STABILITY AND TRANSPORTATION TERMS

Store test:

- in the manufacturer's packaging;
- in a dark, dry place with the temperature from +4 °C to +42 °C;
- in a place protected from mechanical actions (friction, pressure, strokes);
- keep away from the ammonia vapor, moisture and direct sunlight.

When stored at this temperature the device is stable for 24 months. Expiration date is indicated on the label of the package of the test. Transport by any kind of transportation with the temperature from -50 °C to +60 °C, sealed. The transport period must not exceed 1 month.

12. WARRANTY

- The Manufacturer shall remedy all defects discovered in any test (the "Defective Product") that result from unsuitable materials or negligent workmanship and which prevent the mechanical functioning or intended use of the Product including, but not limited to, the functions specified in the Manufacturer's specifications for the Product.
- Any warranty will, however, be deemed as void if fault is found to have been caused by maltreatment, misuse, accidental damage, incorrect storage or use of the product for operations outside their specified limitation or outside their specifications, contrary to the instructions given in the instruction manual.
- The period of this warranty is 24 months from the date of manufacture.

13. ORDERING INFORMATION

- AMA-Med Oy, Sammonkatu 12, 50130, Mikkeli, Finland, tel: +358 45 164 44 04, e-mail: expert@amarut.com
- Manufacturer: Association of Medicine and Analytics Company Limited, 17 line of Vasilevsky Island, 4-6, 199034, Saint-Petersburg, Russia, tel: (007) 812 321-7501, e-mail: ama@sp.ru; www.amamed.eu

14. REFERENCES

- Mobley H.L.T., Mendz G.L., Hazell S.L., editors. *Helicobacter pylori: Physiology and Genetics*. Washington (DC). ASM Press. 2001.
- Ji-Hyun Seo, Ji-Sook Park, Kwang-Ho Rhee, and Hee-Shang Youn. Diagnosis of *Helicobacter pylori* Infection in Children and Adolescents in Korea. *Pediatr Gastroenterol Hepatol Nutr*. 2018 Oct; 21(4): 219–233.
- Lash J.G., Genta R.M. Adherence to the Sydney System guidelines increases the detection of *Helicobacter gastritis* and intestinal metaplasia in 400 738 sets of gastric biopsies. *Alimentary Pharmacology & Therapeutics*. 2013. Aug; 38(4): 424-431.
- Mollenkopf C, Steinger H, Weineck G, Meyer M. Gastritis: immunohistochemical detection of specific and nonspecific immune response to *Helicobacter pylori*. *Z Gastroenterol*. 1990 Jul; 28(7): 327-34.
- Federal State Budgetary Institution "Federal Scientific and Clinical Center for Specialized Types of Medical Aid and Medical Technologies of the Federal Medical and Biological Agency", Clinical Product Evaluation Protocol, Russia, 2019.
- Saint-Petersburg Pasteur Institute of Epidemiology and Microbiology, St. Petersburg, Russia. Study Protocol No. 6/2, 2017.
- Lash J.G., Genta R.M. Adherence to the Sydney System guidelines increases the detection of *Helicobacter gastritis* and intestinal metaplasia in 400 738 sets of gastric biopsies. *Alimentary Pharmacology & Therapeutics*: 2013 Aug; 38(4): 424-431.
- Bode G, Mauch F, Mallertheiner P. The coccoid forms of *Helicobacter pylori* – criteria for their viability. *Epidemiol. Infect.* 1993. 111: 483-490.
- Cellini L, Allocati N, Angellitti D, et al. Coccoid *Helicobacter pylori* not culturable in vitro reverts in mice. *Microbiol Immunol*. 1994. 38: 834-850.
- Takahiro Uotani, David Y. Graham. Diagnosis of *Helicobacter pylori* using the rapid urease test. *Ann Transl Med*. 2015. Jan; 3(1): 9.

EXPLANATION OF THE SYMBOLS, USED IN LABELS

	Manufacturer		Contents sufficient for n tests
	Authorized Representative in Europe		Do not reuse
	Use by		Consult Instructions for use
	Batch code		In vitro diagnostic medical device
	Date of manufacture		Temperature limitation
	European Conformity		Catalogue number

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